BREAST ULTRASOUND IMAGE ACQUISITION PROTOCOL

ENGLISH

Pag. 2

ITALIANO

Pag. 3

BREAST ULTRASOUND IMAGE ACQUISITION PROTOCOL

The emphasis of the study is on the morphometry and echostructure of the breast masses; therefore the image to be acquired is an ultrasound study using a breast ultrasound probe of a "breast mass", defined as a persistent neoformation inconsistent with the physiological breast tissue activity, affecting the breast. The desmoplastic reaction, if present, must be considered as belonging to the lesion. The study lasts about 10 minutes.

Frame N. 1: Display on the longitudinal plane

Correct positioning of the probe on the patient's breast is essential for processing the reflected image of the breast mass with TRACE4BUS.

Please manually move the probe by checking the reflected image of the breast mass on the ultrasound screen.

The patient's lesion should appear in the screen image possibly centered in the reflex image and acquired in the longitudinal section. At this point, make the first still image (frame N. 1) without annotations inside the reflected image and save it in DICOM format.

Please use the still images below as a reference to determine if the lesion is correctly positioned in the still image and free of marks.

The positioning of the mass and the absence of annotations are important for a correct analysis of frame N.1

If the lesion is not positioned correctly in the still image, adjust the positioning of the breast probe by checking the image of the lesion on the screen and re-perform the freeze image. Also make sure that the patient is lying down correctly. Continue repositioning the probe and scanning until the breast lesion is properly centered in the screen image.

Most modern ultrasound systems provide breast probes and automated adjustment procedures for rate adjustments in ultrasound procedures in the breast region. Follow manufacturers recommendations for breast probe selection and adjustment procedures provided by manufacturers.

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