Trace4BUS

PRODUCT SHEET

LANGUAGES

p. 3 ENGLISH

p. 9 ITALIANO

THE PROBLEM

Breast cancer (BC) is one the most frequent malignant neoplasms in women, with 2.3 million women newly diagnosed with this disease and 685 000 deaths worldwide in 2020 according to the World Health Organization.

Ultrasound imaging is a key tool in BC detection and diagnosis. However, benign breast abnormalities can mimic malignancies so that percutaneous needle biopsy is considered the best method for BC confirmation or exclusion to prevent unnecessary surgery, as recommended by the European Society of Breast Cancer Specialists (EUSOMA) and the European Commission Initiative on Breast Cancer.

As a consequence, to reliably exclude malignancy, breast radiologists perform biopsy not only of clearly malignant but also of probably benign lesions, unless they think that a given lesion in a given patient, also considering the patient-specific risk profile, has an extremely low probability of being malignant and that a six-month delayed diagnosis would not impact on patient's outcome.

This clinical behavior results in a huge number of needle biopsy performed every year for breast lesions that are finally revealed to be benign. Considering a rate of benign breast findings on the total number of breast needle biopsies being at least 60% and that at least 70% of them are performed under ultrasound guidance, *every year about one million ultrasound-guided breast biopsy of benign lesions is performed in the world, increasing the psychologic and economic cost of breast care.*


OUR SOLUTION: TRACE4BUS®

The barest masses appear to have characteristics at ultrasonography (US) that are helpful in the diagnosis to expert operators. However, BI-RADS classification guidelines developed so far are helpful in very low risk group discrimination (0% likelihood of malignancy, BI-RADS 2), while showing limits for masses assigned with an intermediate category (BI-RADS 3-4) and complex in some cases for masses assigned with BI-RADS 5, leading to high ratio of false positives.

TRACE4BUS is our proprietary licensed software that uses machine learning to deeply analyze US images of suspicious breast masses of women at risks to provide the breast mass BI-RADS category and malignancy likelihood.

DESCRIPTION

TRACE4BUS provides the diagnostic category of breast suspicious masses according to the BI-RADS classification by a semi-automatic reading of the mass ultrasound images.

INDICATION OF USE of TRACE4BUS®

TRACE4BUS is indicated for subjects with breast masses detected at ultrasound imaging and considered suspicious for breast cancer by physicians with experience in breast ultrasound. Physicians with experience in breast ultrasound can use TRACE4BUS as a support for: reporting ultrasound-detected suspicious breast masses; assigning the diagnostic category to such masses according to the BI-RADS classification; and, as a consequence, making an informed decision to send or not to send such masses to ultrasound-guided needle biopsy. Physicians with experience in breast ultrasound can use TRACE4BUS as a support for diagnosing and managing patients with ultrasound-detected suspicious breast masses. The TRACE4BUS software provides an assessment of suspicious breast masses at the date of the ultrasound examination in which such masses are detected.

It should be noted that DeepTrace Technologies S.R.L. considers TRACE4BUS as a support to physicians with experience in breast ultrasound for reporting, diagnosis and decision-making of suspicious breast masses. However, these medical professionals retain the ultimate decision-making responsibility.

HOW TRACE4BUS® PREDICTS THE RISK OF BREAST CANCER

TRACE4BUS is a standalone software based on US image analytics in the breast region, data mining and machine learning classifiers. It is based on the identification of US image features specifically associated to abnormal morphometric and echogenic pattern within and near a breast cancer mass caused by the breast cancer disease.

With TRACE4BUS, morphometric and echogenic features of a suspicious breast mass are extracted from her US breast image (obtained by breast probe without contrast medium), then automatically analyzed by algorithms of support vector machines, a family of machine learning classifiers that, once properly trained, allow automatically distinguishing between three classes of breast masses of interest, in this case:

-BI-RADS 3 vs

-BI-RADS 4 vs

-BI-RADS 5.

TRACE4BUS uses the features of the morphometric and echogenic pattern of the breast masses of the classes of risk described above, masses of which the histological diagnosis is known, to train its learning machines to distinguish the images of the three classes. TRACE4BUS training is supervised by the histological diagnosis of the masses (malignant disease, benign disease) obtained following biopsy or definitive surgery.

When the US image of the breast mass of a subject at risk for breast cancer is provided as input to TRACE4BUS, the software uses its learning machines to classify that breast mass into one of three classes (<2% likelihood of malignancy vs 2%- 95% likelihood of malignancy vs > 95% likelihood of malignancy).

PERFORMANCES

Performance of TRACE4BUS was tested on a set of US breast images obtained by US systems from different vendors from patients of the SENORETRO study¹ with suspicious breast masses.

In the SENORETRO study, TRACE4BUS was used to predict the risk of malignancy of the breast masses based on the US appearance and was proven to be effective in predicting BI-RADS 3, 4, and 5 classes and in providing an over 15% reduction of the biopsy rate of lesions finally revealed as benign, while still warranting very high sensitivity (>94%) [2].

SENSITIVITY > 94%*

FALSE POSITIVES REDUCTION > 15%

* performance obtained using mono-vendor US systems are:
SENSITIVITY > 98 %, FALSE POSITIVES REDUCTION > 18%

Example of US breast image of a breast mass for a patient at risk of breast cancer, with highlighted the segmented region of interest of the breast mass.

[1] Protocollo SENORETRO IRCCS Policlinico San Donato (San Donato, Italia) v01/17 alla data del 16 December 2021.

[2] Interlenghi M, et al., A Machine Learning Ensemble Based on Radiomics to Predict BI-RADS Category and Reduce the Biopsy Rate of Ultrasound-Detected Suspicious Breast Masses. *Diagnostics (Basel)*. 2022 Jan 13;12(1):187. doi: 10.3390/diagnostics12010187.

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